

i-LIGHT^N®

Light Curing Nano Hybrid Composite



Based on:
UDMA Bis-GMA,
EBPADMA TEGDMA



Zirconium and silica/
barium glass filler:
main glass cluster size
0.25-1.5 µm



Inorganic filler
loading by weight
77,63 %

CLASS IV RESTORATION USING i-FLOW AND i-LIGHT

BEFORE



AFTER



UAB "Medicinos linija"
Address: Aviacijos g. 28, Siauliai
LT-77103, Lietuva
Email: info@i-dental.lt
Website: www.i-dental.lt

i-dental
Reliable solutions for everyday dentistry