

Opalescence® **BOOST**

No light needed

Opalescence Boost in-office whitening is chemically activated, so it doesn't require a light or laser to whiten smiles. It safely provides excellent whitening results in about an hour.

- High water content to prevent shade relapse and dehydration
- Two to three 20-minute applications
- Syringe-to-syringe mixing ensures maximum strength
- Distinct red color aids proper placement and ensures complete removal
- 40% hydrogen peroxide
- Contains PF (potassium nitrate and fluoride)
- No refrigeration required for storage*



ULTRADENT
PRODUCTS, INC.

800.552.5512 ULTRADENT.COM



Before



After

Opalescence Boost in-office whitening delivers visible results in two to three 20-minute applications. Its convenient syringe-to-syringe chemical activation process ensures maximum strength and saves patients from the discomfort of lights or lasers.¹



4750-US Opalescence Boost PF Intro Kit

4 x 1.2 ml Opalescence Boost/Activator syringes
2 x 1.2 ml OpalDam Green syringes
2 x IsoBlocks
2 x Shade guide cards
2 x Ultradent Luer Vacuum Adapters
2 x SST tips
10 x Micro 20 ga tips
10 x Micro 20 ga FX tips



4751-US Opalescence Boost PF Patient Kit

2 x 1.2 ml Opalescence Boost/Activator syringes
1 x 1.2 ml OpalDam Green syringe
1 x IsoBlock
1 x Shade guide card
5 x Micro 20 ga tips
5 x Micro 20 ga FX tips



4754-US Opalescence Boost PF Econo Refill

20 x 1.2 ml Opalescence Boost/Activator syringes



1. Buchalla W, Attin T. External bleaching therapy with activation by heat, light or laser—a systematic review. Dent Mater. 2007 May; 23(5):586-96. *Mixed whitening gel must be refrigerated.

